

Operators: <u>CK, JB, JBD</u>		Test Protocol Subject-Nr.: <u>1208596</u>		Date: <u>27/11/18</u> Time: <u>09:00</u>	
<b>GENERAL PREAMBLE</b>					
- Rooms booked		<input checked="" type="checkbox"/>	- Talkback active (speaker!)		<input checked="" type="checkbox"/>
- Consent forms printed, pen available		<input checked="" type="checkbox"/>	- Fresh earbuds installed with Subject		<input checked="" type="checkbox"/>
- Randomized parameters prepared		<input checked="" type="checkbox"/>	- Chair in correct position		<input checked="" type="checkbox"/>
- Audiometry gear switched On (AMP!)		<input checked="" type="checkbox"/>	- Television running		<input checked="" type="checkbox"/>
- Audiometry run script prepared / comm.		<input checked="" type="checkbox"/>	- General procedure explained to Subject		<input checked="" type="checkbox"/>
- Audio recording started		<input checked="" type="checkbox"/>	- Consent forms signed		<input checked="" type="checkbox"/>
- Affinity in match mode (left/right, levels)		<input checked="" type="checkbox"/>	- Subject and headphones places correctly		<input checked="" type="checkbox"/>
<b>AUDIOMETRY</b>					
- Procedure explained		<input checked="" type="checkbox"/>	- Audiometry left ear		<input checked="" type="checkbox"/>
- Correct Frequencies		<input checked="" type="checkbox"/>	- Audiometry right ear		<input checked="" type="checkbox"/>
			Normal hearing		<input checked="" type="checkbox"/> Yes <input type="checkbox"/> No
			- Familiarisation		<input checked="" type="checkbox"/>
<b>LOUDNESS MATCHING</b>					
- Procedure explained		<input checked="" type="checkbox"/>	- Conductors placed		<input checked="" type="checkbox"/>
			Level for HINT: <u>51</u> dB		
First run:		Second run:		Third run:	
1. match <u>49</u> dB		1. match <u>50</u> dB		1. match <u>50</u> dB	
2. match <u>48</u> dB		2. match <u>52</u> dB		2. match <u>52</u> dB	
3. match <u>51</u> dB		3. match _____ dB		3. match _____ dB	
4. match _____ dB		4. match _____ dB		4. match _____ dB	
<b>HEARING IN NOISE</b>					
- Affinity in HINT mode (Check AC level)		<input checked="" type="checkbox"/>	Familiarisation (Half Air, half bone)		
- BC level set according to matching result		<input checked="" type="checkbox"/>	List number: _____ % <u>2</u> %		
- Conductors placed correctly on Subject		<input checked="" type="checkbox"/>	Threshold <u>12</u> dB		
Test run 1			Test run 2		
Presentation mode		AC <input type="checkbox"/> BC <input checked="" type="checkbox"/>	Presentation mode		AC <input checked="" type="checkbox"/> BC <input type="checkbox"/>
List number:		% <u>5</u> %	List number:		% <u>8</u> %
Threshold		<u>11</u> dB	Threshold		<u>11</u> dB
Test run 3			Test run 4		
Presentation mode		AC <input type="checkbox"/> BC <input checked="" type="checkbox"/>	Presentation mode		AC <input checked="" type="checkbox"/> BC <input type="checkbox"/>
List number:		% <u>3</u> %	List number:		% <u>9</u> %
Threshold		<u>84</u> dB	Threshold		<u>11</u> dB
What transducer was perceived to be more challenging for HINT?				AC	BC <input checked="" type="checkbox"/> n/a
Comments: <u>Subject has hard time differentiating BC and AC.</u>					
<b>DEBRIEFING</b>					
- Debrief Subject on the test		<input checked="" type="checkbox"/>	- Organize paperwork		<input checked="" type="checkbox"/>
- Provide coffee, snacks and candy		<input checked="" type="checkbox"/>	- Backup test data/recordings on USB stick		<input checked="" type="checkbox"/>
- Remove dirty Earbuds		<input checked="" type="checkbox"/>	- Test protocol completely filled out		<input checked="" type="checkbox"/>
- Want to know end results?		<u>YES</u>			